



H3ABioNet

Pan African Bioinformatics Network for H3Africa

16SrRNA Intermediate Bioinformatics Online Course: Int_BT

Module 3:

Sample collection, extraction and library preparation for 16S NGS analyses



H3ABioNet

Pan African Bioinformatics Network for H3Africa



16SrRNA Intermediate Bioinformatics Online Course:

Int_BT_2019

Shantelle Claassen-Weitz

Learning Outcomes

Describing different aspects of planning for 16S experiments; for example study design, DNA extraction methods and laboratory workflows.

- Understand what is meant by the term “16S rRNA gene” and why we are interested in this gene.
- Understand the concepts behind Sanger sequencing and high-throughput sequencing
- Know why it is important to clearly plan an experiment and how different components to the study and experimental design may influence data generated and downstream analyses.
- Know that techniques other than 16S rRNA sequencing are also available to study microbial profiles.



H3ABioNet

Pan African Bioinformatics Network for H3Africa

16SrRNA Intermediate Bioinformatics Online Course: Int_BT_2019

Module 3:

Sample collection, extraction and library prep for 16S NGS analyses

Part 3.1

16S rRNA high throughput sequencing: how it works



H3ABioNet

Pan African Bioinformatics Network for H3Africa



16SrRNA Intermediate Bioinformatics Online Course:

Int_BT_2019

Shantelle Claassen-Weitz


16S rRNA high throughput sequencing: how it works




mid-1600s:
First microbes described




1800s – Present:
Culture, staining, and microscopy used to study microbes that can be cultured



circa 1600:
Microscope invented



1800s:
Connection made between microbes and disease



1990s:
DNA sequencing becomes available, allowing study of microbes that cannot be cultured

Advantages: Study viable organisms
Disadvantages: Thought that less than 1% of all bacterial species are cultivatable

16S rRNA high throughput sequencing: how it works

Culture-independent molecular methods:

Marker gene-dependent



Fingerprinting techniques

Sanger sequencing

High-throughput sequencing



16S rRNA gene

All organisms need ribosomes to make protein

Any gene that makes up ribosome may be a good marker gene



Ribosomal RNA (rRNA) never gets transcribed to protein

Gene focussed on: 16SrRNA gene (~1500bp)

Marker gene-independent





High-throughput sequencing (shotgun sequencing)



https://www.youtube.com/watch?v=8Aa_mnyXm70

16SrRNA Intermediate Bioinformatics Online Course:
Int_BT_2019
Shantelle Claassen-Weitz

16S rRNA high throughput sequencing: how it works




The entire gene is too long to sequence without errors for most DNA sequencers

Find a portions with enough variability to distinguish species



https://www.youtube.com/watch?v=8Aa_mnyXm70

16SrRNA Intermediate Bioinformatics Online Course:
Int_BT_2019
Shantelle Claassen-Weitz


16S rRNA high throughput sequencing: how it works




16S ribosomal RNA gene



DNA amplicons




Sanger method only sequences a **single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works

Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works

Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works

Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works

Sanger method sequences **a single DNA fragment at a time**




Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works

Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works


Sanger method sequences **a single DNA fragment at a time**



Sanger method sequences **a single DNA fragment at a time**




Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works

Sanger method sequences **a single DNA fragment at a time**




16S rRNA high throughput sequencing: how it works






16S rRNA high throughput sequencing: how it works


DNA from a number of samples



Undergo amplification using a multiplex approach



16S rRNA high throughput sequencing: how it works



In summary:

- 16S rRNA gene is by far the most common housekeeping gene targeted to study bacterial phylogeny and classification.
- High throughput sequencing targeting the 16S rRNA gene allows for large quantities of DNA to be sequenced much more quickly and cheaply compared to Sanger sequencing.
- High throughput sequencing of the 16S rRNA gene allows for identification as well as relative abundance quantification of all bacteria present in a sample.
- 16S rRNA high-throughput sequencing allows for processing multiple samples together via the use of barcoded (indexed) primers.